



GI RESTORATIVE

SECURE GI FILL XTRA 9

THE DURABLE AESTHETIC GI RESTORATIVE

Secure GI FILL XTRA 9 packs

Extraa Strength

Extraaaa Durability

Extraaaaa Translucency

Extraaaaaa Aesthetics and

Extraaaaaaa Savings in a single pack.

Incorporating a fine balance of strength and aesthetics, this GI filling material works perfectly for anterior as well as posterior restorations. Thus, helping reduce inventory and cost.



wizdent





Secure GI FILL XTRA 9 specially ground powder particles make it super easy to mix and give the mix a nice packable, non-sticky consistency for effortless handling.

Best in class **Compressive strength (258 MPa)** helps deliver durable long-lasting restorations.

Added Calcium Oxide imparts **extra translucency** resulting in superior aesthetics



Unique bridge shade concept incorporates **2 shades (A2-A3) in a single bottle**, thereby saving cost

Excellent biocompatibility, Chemical bonding and low solubility ensure **excellent marginal integrity**

Quick mixing and **easy handling** make it an ideal material for patients of all ages.

It's liberal fluoride release feature helps prevents plaque build-up and reduces secondary caries.



SPECIALLY INDICATED FOR

Class I *

Class III

Core Build Up



Class II *

Class V

*Non stress bearing

USAGE

- ⌚ Mixing Time 30 secs
- ⌚ Working Time** 1 min 30 secs (minimum)
- ⌚ Setting Time*** 3 mins (maximum)

Working Time is at $23 \pm 1^{\circ}\text{C}$, *setting time is at $37 \pm 1^{\circ}\text{C}$

ORDERING INFORMATION

- 📦 15 gm powder
- 📦 7 ml liquid
- 📦 1 shade – A2 ••• A3

[Click here](#) to find more on what other dentists are talking about Wizdent Products

📞 7506698711

✉ wizdent@pidilite.com

📘 fb.com/fineartofdentistry

🌐 www.wizdent.com