



GLASS IONOMER **SECURE CORE Z**

THE TRUE SELF ADHESIVE BULK FILL

For years, practitioners have traded off the biocompatibility, anti-cariogenic properties of GLC's and marginal integrity for the strength of composites filling materials. Secure Core Z combines the benefits of both into a single product.

Reinforced with Zirconia fillers Secure Core Z delivers "Rock Solid" restorations that bond chemically to the tooth.





This ingenious Zirconia infused self-adhesive restorative has been designed specifically for durable, single step bulk restoration - Class I,

Outstanding strength - compressive strength similar to that of composite filling materials 346 MPa

Molecular bonding, high biocompatibility and **fluoride release** alike glass ionomers ensures better marginal integrity

Packable consistency allows easy **one step bulk placement** eliminating the need for incremental placement

Self-cure feature minimizes shrinkage and virtually eliminates post-operative hypersensitivity

Unique A2-A3 **bridge shade** helps minimize inventories

Secure Core Z helps you deliver durable, long lasting posterior restorations with great ease and confidence.



SPECIALLY INDICATED FOR

- Bulk restoration of Class I and Class II cavities in permanent & deciduous teeth
- Core build-up restorative under crowns
- Amalgam replacement
- On the root surfaces for locating overdentures

USAGE

- ⌚ Mixing Time 30 secs
- ⌚ Working Time** 1 min 30 secs
- ⌚ Setting Time*** 3 mins 30 secs

Working Time is at $23 \pm 1^{\circ}\text{C}$, *setting time is at $37 \pm 1^{\circ}\text{C}$

ORDERING INFORMATION

- 🧴 15 gm powder
- 💧 10 ml liquid
- 1 shade – A2 ••• A3

[Click here](#) to find more on what other dentists are talking about Wizdent Products

📞 7506698711

🌐 fb.com/fineartofdentistry

✉ wizdent@pidilite.com

🌐 www.wizdent.com