



LUTING CEMENT

# WONDER LUTE UNIVERSAL

## THE ONE FOR ALL SELF ADHESIVE RESIN CEMENT

The recent developments in crown & bridge materials and the available cementation options have made routine cementation procedures complex and confusing. Which cement works best with which substrate? What steps need to be followed? Is surface treatment required? Which primer should be used? Questions like these are baffling clinician's the world over.

Wonder Lute Universal helps tackle this problem with it's unique, super easy to use formulation.



**wizdent**





An affordable everyday luting cement, Wonder Lute Universal requires no separate etching /priming. This highly aesthetic, **“One For All”** resin cement is the ideal choice for all kinds of permanent indirect restorations – metal, zirconia, all ceramic. Lithium disilicate & post cementation.

**Super Easy application** - Just mix, apply and cure!

**Highly Affordable** –saves time & cost

High Bond strength ensures **superior retention**

Virtually **insoluble** with **no post-operative sensitivity**

**Easy clean up** with tack cure feature

Low film thickness of 12µm enables **accurate fit**


**Dual cure** with command set feature





Wonder Lute Universal is available in an aesthetic translucent shade and is also radiopaque for easy identification on radiographs.



## SPECIALLY INDICATED FOR




All post  
cementation


Crown / Bridge / Inlay  
Cementation

## ORDERING INFORMATION

-  5ml Dual Barrelled Syringe x 1
-  1 Shade – Translucent

[Click here](#) to find more on what other dentists are talking about Wizdent Products

 7506698711

 fb.com/fineartofdentistry

 wizdent@pidilite.com

 www.wizdent.com