



LUTING CEMENT **SECURE LUTE 1**

THE RELIABLE GI LUTING CEMENT

Frequent complaints about crown dislodgement are worrisome for both the dentist as well as the patient. Secure Lute 1 is a reliable glass ionomer luting cement inspired by six decades of trusted Pidilite legacy in adhesive technologies.

Designed to overcome challenges like crown dislodgment and secondary marginal caries, Secure Lute 1 promises to make everyday cementation quick and stress-free.





A highly biocompatible luting cement, Secure Lute 1 offers superior compressive strength (128 MPa) and ensures durable, long - lasting cementation with a variety of substrates.

• **Pre-dosed scoop** calibrated to Indian teeth sizes that minimizes material wastage and ensures a consistent mix every time

• **Low Film Thickness** (17 microns) ensures better fitting prothesis



• Extremely **Low Solubility** for Long-lasting cementations

• Unique **Quick Set** formulation that minimizes chair-time

• **Extra liquid** also helps extend the spectrum of usage including cavity conditioning and orthodontic procedures.



SPECIALLY INDICATED FOR



Crown / Bridge
cementation – Zirconia,
PFM and Metal



Cementation of
Inlays / Onlays,
Orthodontic bands
and appliances

USAGE

⌚ Mixing time 20 secs

⌚ Working time** 1 min

⌚ Setting time*** 2 mins

Working Time is at $23 \pm 1^{\circ}\text{C}$, *setting time is at $37 \pm 1^{\circ}\text{C}$

ORDERING INFORMATION

📦 35 g powder

💧 21 ml liquid

[Click here](#) to find more on what other dentists are
talking about Wizdent Products

📞 7506698711

🌐 fb.com/fineartofdentistry

✉ wizdent@pidilite.com

🌐 www.wizdent.com