



MASTER DESIGN

Deliver high end aesthetic restorations swiftly & with virtually no instrumentation

Aesthetic composites restorations generally take time and require extensive armamentarium which can be challenging in situations like the ongoing COVID 19 pandemic.

Master Design, a speciality nano-restorative enables quick, free-hand direct restoration of teeth with virtually zero instrumentation. It's stackable, non-slumping consistency and controlled fluidity feature help minimize patient chair-time and the risk cross-contamination.

- └ Uniquely engineered consistency for **Instrument Free Adaptive Dentistry**
- └ **High compressive strength (340 MPa)**, ensures long-lasting posterior restorations
- └ Virtually **void free restorations** for improved strength and durability
- └ Experience the comfort of **direct material stacking and anatomy recreation** with the composite syringe itself
- └ Specifically suited for quick construction **proximal walls, cervical hip, palatal shelves, deep margin acquisition**
- └ Unparalleled **ease and comfort**





SPECIALLY INDICATED FOR

- Direct Class I, II, III, IV, and V restorations
- Reconstructing the proximal box and gingival seat of Class II case
- Palatal wall reconstruction in Class III & Class IV
- The characterization in anterior aesthetic procedures, e.g. Mamelon creation, Recreation of incisal translucency
- Diastema closure and Occlusal anatomy recreation
- Sealing hypersensitive areas and blocking undercuts
- Splinting mobile teeth

ORDERING INFORMATION



3 g Syringe x 1



Applicator Tips x 5



2 Shades – A1B A2B A3B A3.5B B1B
B2B LD MD DD TE

DENTAL MATERIALS FOR SAFER PRACTICES

CORE BUILD-UP



GI LUTING CEMENT



BULK FILL COMPOSITE



LIGHT CURE RESTORATIVE GIC



SPECIALITY COMPOSITE

