

QUICK SET!



SECURE LUTE 1

A reliable GI Luting Cement that helps minimize cementation time

Everyday Type 1 Luting cements have a self-cure setting time of more than 2 minutes which is worth considering amidst the ongoing Covid-19 pandemic. Quick to mix and apply, Secure Lute 1 promises to make everyday cementation **fast and stress-free**. Its Unique Quick Set formulation helps reduce cementation time **by up-to 50%**.



Secure Lute 1 is a reliable GI Luting Cement inspired by six decades of trusted Pidilite legacy in adhesive technologies.

- Low Film Thickness (17 microns) ensures better fitting prosthesis.
- Superior compressive strength (128 MPa) for durable, long - lasting cementation with a variety of substrates.
- Extremely Low Solubility for Long-lasting cementations.
- Pre-dosed scoop calibrated to Indian teeth sizes minimizes material wastage and ensures a consistent mix every time.
- Extra liquid ensures you don't have to purchase a separate bottle of liquid.





SPECIALLY INDICATED FOR



Crown / Bridge
cementation – Zirconia,
PFM and Metal



Cementation of
Inlays / Onlays, Orthodontic
bands and appliances

USAGE



Mixing time 20 secs



Working time^{**} 1 min



Setting time^{***} 2 mins

^{**}Working Time is at $23 \pm 1^\circ\text{C}$, ^{***}setting time is at $37 \pm 1^\circ\text{C}$

ORDERING INFORMATION



35 g powder



21 ml liquid

DENTAL MATERIALS FOR SAFER PRACTICES

CORE BUILD-UP



GI LUTING CEMENT



BULK FILL COMPOSITE



LIGHT CURE RESTORATIVE GIC



SPECIALITY COMPOSITE

