

A close-up photograph of two hands against a dark background. The hands are holding several thin, white threads that hang down. The lighting is dramatic, highlighting the texture of the skin and the threads.

SECURE GI FILL LC

A durable resin reinforced GI restorative that sets on your command!

Conventional Restorative GIC's take minimum 3 minutes for self-curing. This precious time can be reduced to seconds while treating patients during this pandemic, when our priority is to reduce chair-side time, thereby reducing risk exposure.



SECURE GI FILL LC

SECURE GI FILL LC's unique command cure feature cures the GI restoration in just 20 secs (upto 3mm). Its non-sticky, packable consistency makes it quick and easy to place.

- Its specially incorporated dual cure feature allows adequate working time to shape & mold the restoration – reducing the need for finishing & minimizing aerosol production.
- It also allows you to finish and polish the restoration in the same sitting – no repeat visits.
- Non sticky, packable consistency minimizes need for instrumentation





SPECIALLY INDICATED FOR

Class I *

Class III

Core Build Up







Class II *

Class V

*Non stress bearing


USAGE

-  Mixing Time 30 secs
-  Working Time** 2 min (minimum)
-  Self Cure Setting Time*** 6 mins
-  Light Cure 20 secs (minimum)

Working Time is at $23 \pm 1^\circ\text{C}$, *Self Cure Setting is at $37 \pm 1^\circ\text{C}$

ORDERING INFORMATION

 15 gm powder

 6 ml liquid

•• 1 shade – A2 ↔ A3

Dental Materials For Safer Practices

CORE BUILD-UP



LUTING CEMENT



BULK FILL COMPOSITE



LIGHT CURE RESTORATIVE GIC



SPECIALITY COMPOSITE

