



COMPLEX

SIMPLE

WONDER BULK FLOW

A Self Adaptive Bulk Fill that makes complex posterior restorations Quicker & Simpler

During these pandemic times, complex posterior restorations need smart composites that help reduce armamentarium and chair-side time to avoid unnecessary exposure. Ease of application without compromising on clinical efficacy is vital ...



Wonder Bulk Flow is an advanced self-levelling adaptive bulk fill composite designed to deliver one step curing depth of up-to 4mm. This makes Wonder Bulk Flow an ideal choice for posterior fillings and tunnel restorations.

- └ Convenient and easy adaptation with self-levelling properties
- └ Quick and stress-free bulk placement with one-time cure of up-to 4mm
- └ Best in class filler load of 73% ensures high compressive strength with no capping layer required



Equipped with special **symmetric refractive fillers**, no need to compromise on shade matching while using a bulk composite. Moreover, these fillers intelligently adapt to the optical density of surrounding tooth structure, hence achieving an appropriate shade match with a single universal syringe.

SPECIALLY INDICATED FOR



UP TO 4 mm
SINGLE STEP CURE

Bulk Posterior
Restorations



Tunnel
Restoration



Class V
Abfractions

ORDERING INFORMATION

— 2g Syringe x 1

— Applicator Tips x 5

•• 1 Shade – Universal

DENTAL MATERIALS FOR SAFER PRACTICES

CORE BUILD-UP



GI LUTING CEMENT



BULK FILL COMPOSITE



LIGHT CURE RESTORATIVE GIC



SPECIALITY COMPOSITE

